

## Experiment 2

---

# Geometrical Optics

---

### 1 Introduction

In this experiment, we will continue to explore geometrical optics by studying the optics of simple curved mirrors and lenses.

### 2 Background - see Pedrotti<sup>3</sup>, Sections 2-6 to 2-9

When studying the geometrical optics of mirrors or lenses one considers the following three quantities: object distance  $s_o$ , image distance  $s_i$ , and focal length  $f$ . These quantities are related by the equation

$$\frac{1}{s_o} + \frac{1}{s_i} = \frac{1}{f} \quad (1)$$

There is a convention to be followed in the definition of these quantities. For lenses, a converging lens (convex) has  $f > 0$  while a diverging lens (concave) has  $f < 0$ . For mirrors,  $f > 0$  for concave mirrors, and  $f < 0$  for convex mirrors. Also by convention, we place the object to the left of the lens, with  $s_o > 0$ . If  $s_i > 0$ , it is on the right of the lens and is a *real* image. If  $s_i < 0$  it is to the left of the lens (same side as object) and is a *virtual* image. One can consider the mirror as a folded over version of the lens:  $s_o$  is positive and on the left, but now a  $s_i > 0$  is on the left (the opposite of the lens) and  $s_i < 0$  is on the right, behind the mirror, and a virtual image.

The focal length of a spherical mirror is simply  $f = R/2$ , where  $R$  is the radius of the mirror, and the focal length of a thin lens is given by

$$\frac{1}{f} = (n - 1) \left( \frac{1}{R_1} - \frac{1}{R_2} \right), \quad (2)$$

where  $n$  is the index of refraction, and  $R_i$  are the radii of curvature of the two surfaces.

### **3 Experiment**

You are supplied with a concave mirror, and two lenses. Your challenge is to find the focal lengths of the optics as accurately as possible.

A first start is to use an object at infinity (or close to it). The windows in the hallway are a good place to start.

You have available the laser as a source of rays to do ray tracing. Think about what happens to the rays as they pass through a lens, or bounce off a mirror. Also available is a scanning photodiode, which will allow you to use the computer to acquire beam position and sizes.

Using two (or more) optics, build a compound lens device, and measure its properties.