



PHYS 275 – Experiment 8

Standing Waves on a String




Experiment Summary



- Today we will study oscillations in a string under tension
 - The physics is a bit more complicated than last weeks: we will study standing waves
 - The experimental apparatus consists of a string under tension (a weight is applied to it); we will measure the frequency of its oscillations
- Technical topic of the week: weighed averages
 - We will combine measurements using weights, which will allow us to improve the precision of our estimate of the quantity of interest

Standing-waves Theory



- Standing (or stationary) waves are distinguished using their mode number n
 - Node = location where the amplitude is zero
 - Mode number n = number of nodes + 1
 - Wavelength $\lambda = 2L/n$

More Theory

- The speed and frequency of a standing wave are related:

$$- v_s = \lambda f = \sqrt{\frac{T}{\mu}}$$

This is a form of
dispersion relation

- T is the tension of the string, measured in Newton
- μ is the linear mass density of the string (mass/length), measured in kg/m
- λ is the wavelength, $\lambda = \frac{2L}{n}$
- f is the wave frequency, measured in ...
- v_s is the wave speed

$$\Rightarrow \frac{f_n}{n} = \left[\frac{1}{2L\sqrt{\mu}} \right] T^{1/2}$$



Measurements

- We will measure the frequency corresponding to different modes, and combine these measurements to obtain the best estimate of the linear mass density
 - The idea is that when you combine measurements with different precision, you want to give more weight to the more precise ones
- Here are the main formulae
 - Let us assume that we have the measurements $(x_1, \sigma_1), (x_2, \sigma_2), \dots, (x_n, \sigma_n)$
 - We take the weight of each measurement to be $w_i = \frac{1}{\sigma_i^2}$
 - Finally: $\langle x \rangle = \frac{\sum_{i=1}^n w_i x_i}{\sum_{i=1}^n w_i} = \frac{\sum_{i=1}^n x_i / \sigma_i^2}{\sum_{i=1}^n 1 / \sigma_i^2}$; $\sigma_{\langle x \rangle} = \frac{1}{\sqrt{\sum_{i=1}^n w_i}} = \frac{1}{\sqrt{\sum_{i=1}^n 1 / \sigma_i^2}}$

Experimental Apparatus

Variable frequency & tension standing wave setup:

Frequency generator;
Permanent magnet;
string; meter stick;
weight holder; slotted weights;



Some Suggestions




Use this meter to read frequency

Don't peel off this black tape.
Frequency reading from this
panel is not accurate

Be careful. This
resistor can
become very **hot!**

Adjust magnet position and current
amplitude to make sure string
doesn't touch magnet




More Suggestions

- mass density, μ

$$\mu = M/L$$

We have sample string for you to measure mass density. – everyone needs to measure mass of sample string; $L = 2.0000 \pm 0.0003\text{m}$

- Determine the uncertainty of f_n : (1) the difference between the resonant frequency and the frequency at which the rapid increase in amplitude begins



or (2) measure frequency for a few times and calculate standard deviation as uncertainty.

- pay attention to UNITS

(Yet) More Suggestions

- Calculation of weight and tension: need to remember to **plus the mass of the holder**. Alternatively, my suggestion is to bring your mass with mass holder and to measure together by scale
- Vertical mode standing wave: try to decrease the amplitude; but if you still can not get rid of circular mode, please don't worry- you can simply take it as systematic error.
- Calculation of weighted (f_n/n): I will suggest making a table to calculate all the components you need for weighted mean:

n	f_n	f_i	$\sigma(f_n)$	(f_n/n)	$\sigma(f_n/n)$	$w_n=1/\sigma^2$	$[(f_n/n)] \cdot w_n$
1							
...							

$$\sum_n w_n \quad \sum_n [(f_n/n) \cdot w_n]$$



Notes and Reminders

- Submit your Excel spreadsheet on ELMS and turn in your check sheet before leaving the lab
- Complete the final version of your report by 1pm next week
- Finish the homework set in Expert-TA by 2pm next week
- Turn off your equipment and clean up your bench area before leaving the classroom
- Save your data on the local disk frequently!