



Figure 4: Schematic illustration of how to use a spherometer to measure the radius of curvature of a lens. Note,  $h$  measures the difference in the lengths of the center screw and the three feet and  $b$  is the length of the side of the triangle formed by the feet. **DO NOT USES THE SPHEROMETER ON THE MIRROR!**

**DO NOT USES THE SPHEROMETER ON THE MIRROR!**

Use the spherometer to measure the radii of curvature of each surface of each of the lenses. The configuration is sketched in Fig. 4. First, measure  $b$ , the separation of the legs. Next, calibrate the zero point of the dial and screw on a flat surface. Then, determine the radius of curvature  $R$  from a measurement of the height  $h$  by using the equation

$$R = \frac{b^2/3 + h^2}{2h}.$$

$$\frac{1}{f} = (n - 1) \left( \frac{1}{R_1} - \frac{1}{R_2} \right),$$

relates the focal length to the index of refraction of the lens medium,  $n$ , and to the radii of curvature of the front (entrance),  $R_1$ , and back,  $R_2$ , surfaces of the lens. The following sign conventions have generally been adopted for lenses. For a double concave lens  $R_1 < 0$ ,  $R_2 > 0$ , and  $f < 0$ . For a double convex lens  $R_1 > 0$ ,  $R_2 < 0$ , and  $f > 0$ . If one side of the lens is flat this corresponds to  $R = \infty$ .