Lab 4: How large a force is needed to transition from random to directed motion? Random vs. Directed Motion.

Overview

Much of the first month of Physics 131 focused on directed motion. In both the lecture and laboratory portions of the course, we saw that when the forces applied to objects did not vary substantially from moment to moment, the objects accelerated/established motion in one particular direction.

Since then, we have begun to explore situations in which the pushes and pulls from the environment surrounding the object fluctuate rapidly in time—so rapidly that directed motion is no longer observed. We describe the motion of objects experiencing this volatile bombardment as "random" motion. In lab we have been exploring some of the qualitative and quantitative features of this random motion.

While "random" motion and diffusion are important mechanisms, they propagate very slowly over long distances. Therefore, in many biological processes both random and directed motions are utilized. Over the next two weeks we will take a look at situations where *both* random and directed motion can be observed at the same time. We will try to determine the conditions under which one or the other type of motion is dominant.

Investigation

Our overall goal here is to create conditions in which our beads undergo both random and directed motion, and to find a way to characterize that situation.

We've provided you with a wooden/metal "step" that allows you to tilt the microscope to a couple of different angles. **Be very careful when tilting the microscopes!** Tilting the microscope changes the direction of the normal force exerted by the microscope slide on the beads. What

happens to the motion of the beads when the microscope is tilted in this way?

Because of the dangers of fluid flowing onto the microscope objective, we will be using chambered slides for this lab. Be sure that you fill the chamber to the point where fluid covers the entire surface. This way, we can minimize the amount of fluid flow that occurs, and look only at the motion of the beads relative to the fluid.



- 1. Video #1: Record a video of 5-micron silica beads in water with the microscope tilted using the middle step of the block. A video of 2 to 5 seconds length with about 10 frames per second would be sufficient. Can you observe directed motion? If so, take a video for a long enough time interval that directed motion is clearly visible. Measuring the position of beads for at least 5 time points, create a plot that shows the dependence of r² on time interval for the beads. (If you are using automatic tracking (multitracker plugin), then use all produced data for all 'real' bead tracks—not just 5-10 time points.) Is the dependence of r² on time linear, as it was for the random motion when the microscope was not tilted?
- 2. Video #2: Next, take a sample containing 2-micron silica beads. Record a video of the bead's motion over the <u>same</u> time interval. Track the beads and plot the dependence of r² on time interval. Can you see the directed motion at all or are the beads merely displaying random motion?
- 3. Video #3: Now record a video of the 2-micron silica beads' motion over a larger time interval (a minute or two). For the sake of video size, reduce the frame rate to ~1 fps. Track the beads (using 10 to 15 time points) and plot the dependence of r² on time interval. What type of motion is visible now?

Interpretation

If the motion you are observing in some fluid is *entirely* random, how would you expect r^2 to vary as a function of the time interval over which you measure? If the motion you are observing in some fluid is *entirely* directed, how would you expect r^2 to vary as a function of the time interval over which you measure?

Can a log-log plot help you distinguish between random and directed motion? Does the slope of the log-log plot change at all? Does the slope clearly indicate random or directed motion?

Is it possible to have a situation where both types of motion occur? Which one would 'win'? At what time scales? What are the conditions under which the motion of an object appears to be purely random? What are the conditions under which the motion of an object appears to be purely directed?

Does the concentration of beads that you see on the slide when you look at it under the microscope change at all over the course of the experiment? Can you explain this effect using your results?

¹ Make sure to record the ACTUAL frames-per-second rate used by VirtualDub—the ACTUAL rate is not necessarily the same as what you told the program to do. VirtualDub will examine parameters such as ambient light to determine the best frame rate for the video it collects and will then report the average frames-per-second that it has used once the video has been recorded. Once you move on to collecting/capturing other videos, this frame rate is lost, so it must be recorded as soon as a capture is complete.

² A 'real' track exists if the bead being followed: a) does not enter/leave the selected area over the course of the entire video; b) does not merge with or separate from another bead/track over the course of the entire video; c) does not blink in and out of existence; d) does not 'jump' (track number switching beads, causing a major discontinuity in x- and y-pixel locations); and e) exhibits some form of motion (i.e., is not stuck to the slide/chamber).